GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

Rapor Tarihi Report Date	: 25.08.2024	Son Güncelleme No Last Update No	: 02
SSCP Referans Num. SSCP Reference Num.	: TR-MF-000020130-SSCP02	Son Güncelleme Tarihi Last Update Date	: 05.01.2026

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This summary of safety and clinical performance (SSCP) is intended to provide public access to an updated summary of key aspects of the safety and clinical performance of the device

The SSCP is not intended to replace the instructions for use as the main document to ensure the safe use of the device and is not intended to provide diagnosis or treatment recommendations to target users or patients

The following information is intended for users/healthcare professionals. This information is followed by a summary for patients

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Selecting Screws:

- Bone fixation screws, tent screws and membrane fixation pins are used for directed bone regeneration (GBR) purposes. Tent screws prevent collapse of the membrane placed over the bone grafts. Bone fixation mini screws are used to fix the block bone graft to the donor site. In summary, bone fixation mini screws are indicated when block bone is used, and tent screws are indicated when granular bone graft such as powder/chips + membrane is used. Membrane pins (the smallest mini screws) are used to fix the membrane covering the block or granular bones used in the GBR method. Depending on the physician's preference, these screws and pins can be used together or separately in individual cases.
- If intermaxillary fixation is intended, IMF (Intermaxillary Fixation) mini screws are used.
- Anchorage mini screws are used to move the teeth.

Surgical Technique:

Mini screws are grouped below according to their application purposes;

Membrane Fixation Screws have a flexible driver with a flexible tip that carries them to the operation area and is used to fix them to the bone. On a hard surface, this driver is brought on the pin and pressed at an angle of 90 degrees and the pin is removed from the pin slot with a "click" sound. The screw (pin) at the end of the driver is brought to a certain point of the membrane to be fixed and sent into the bone by hand force or by hitting the lower part of the driver with a hammer. After the screw is inserted into the jawbone, the driver is slightly rotated to the right and left and the screw is left on the membrane by stretching the driver tip. After the time required for GBR (3-10 months) is completed, these screws (pins) are usually removed. Since they are made of biocompatible material, some physicians do not remove these pins at all and leave them in place. There is no medical disadvantage in leaving the pins in the mouth.

Tent Effect Mini Screws are designed to protect the space under the membranes with its self-screwing tip, polished neck and wide head without the need for drilling. In this way, it is aimed to maintain the volume of the bone graft by preventing membrane collapse. It is taken from the set by holding it directly with the driver and fixed to the bone by turning clockwise with light hand pressure. It is applied in such a way that the grooves at the end of the tent screw remain inside the bone and the shiny neck and head parts remain outside the bone. There are driver tips that can be used manually and with contra-angle handpiece. If the tip compatible with the contra-angle handpiece is preferred as the driver during placement, the screws are sent to the bone clockwise (forward) using a speed of 30 rpm and a torque of 20 Ncm.

Bone Fixation Mini Screws are used to fix the block bone graft to the donor site by drilling beforehand. Both the block bone graft and the bone in the donor site are drilled simultaneously with the drill suitable

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

for the screw diameter. The bone fixation screw is inserted into the prepared hole by turning clockwise with the help of a driver (screwdriver) and the block bone graft is stabilized in the desired area. Depending on the size of the bone graft, more than one screw can be used to prevent the graft from moving. There are driver tips that can be used manually and with contra-angle handpiece and are the same as those used in tent screws. If the contra-angle handpiece compatible tip is preferred as the driver during placement, the screws are sent clockwise (forward) into the bone using a speed of 30 rpm and a torque of 10 Ncm.

Intermaxillary Fixation (IMF) Screws. All IMF screws are self-drilling. The risk of damaging dental roots is thus reduced as there is no need for pre-drilling. The ideal position is above the teeth roots.

Using the screwdriver, the screw is directly inserted into the bone; pre-drilling is not necessary. The screw is advanced until the concave shoulder of the screw head is 1 mm from the mucosal surface. The screw may be readjusted for proper alignment of the wire holes contained within the screw head.

As an alternative to the above procedure, the mini screw may be also inserted by using the contra-angle handpiece. In case of insertion by using a contra-angle handpiece, a speed of 30 rpm and a torque of 30 Ncm must be set.

The wires are inserted through the holes of the maxillary and opposing mandibular screws. Alternatively, the wires can be wrapped around the wire groove contained within the head of the IMF screw. The number of wires placed, their configuration, and the amount needed to stabilize the occlusion and fracture is determined by the surgeon. After having established the occlusion, the wires are tightened. The IMF screws may be advanced further to increase the tautness of the wires. Finally, the wires are cut and bent under to prevent soft tissue irritation.

Anchorage Orthodontic Mini Screws. The tip of the screwdriver must be exactly perpendicular to the screw. This position guarantees the retention of the mini screw inside the screwdriver and allows the mini screw to be seized. After having found the engagement, the screwdriver is pressed down to promote the retention of the mini screw inside the tip. The mini screw is inserted to the bone in a clockwise direction. The mini screw is screwed until its head meets the cortical bone of the alveolar ridge. Soft tissues perfectly adhere to the head of the mini screw. The mini screws are self-tapping. In order to screw the implant a pressure is exerted on the knob of the screwdriver with the palm of the hand and the instrument is rotated with the fingers.

As an alternative to the above procedure, the mini screw may be also inserted by using the contra-angle handpiece. When insertion using an contra-angle handpiece, the speed and torque should be adjusted according to the screw diameter as shown below.

Diameter	Insertion torque settings	Insertion Speed
Ø1.40 mm	15 Ncm	30 rpm
Ø1.65 mm	20 Ncm	
Ø2.00 mm	30 Ncm	

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/--

The mini screw can be loaded with chains, elastic wires, or springs for the connection to the orthodontic appliance to induce the required tooth movements. The wire, the chain or the coil spring are connected to the mini screw by using the hole through the head or the groove of the mini screw. The application time of the orthodontic traction on the mini screw depends on the clinician's judgement; the mini screw can usually be loaded immediately after insertion or after healing of soft tissues.

Removing The Screws: Bone fixation mini screws, tent-effect mini screws, IMF mini screws and orthodontic mini screws are removed by turning them counterclockwise with the same driver 1-12 months after insertion. During removal, if the contra-angle handpiece compatible tip is preferred as the driver, the screws are removed counterclockwise (reverse) using a speed of 30 rpm and a torque of 15-30 Ncm. Membrane pins are usually removed after the 3-10 months period required for GBR is completed. No special surgical equipment is required for removal. Periosteal elevator or probe tip is used to gently dislodge the pin. It is removed by holding it with tissue clamp or hemostat. Since they are made of biocompatible material, some physicians never remove these pins and leave them in place. There is no medical disadvantage in leaving the pins in the mouth.

Effect of the Device:

➤ Orthodontic Mini Screw Effect Type:

The mini screw creates a stable anchorage point against orthodontic forces by mechanically attaching to the bone tissue. It exerts its effect through mechanical anchorage and load transfer.

Relation to Medical Device Definition: Orthodontic mini screws are implantable passive medical devices that have no pharmacological effect and operate on the principle of mechanical action.

➤ Intermaxillary Fixation (IMF) Screw Effect Type:

The screw is mechanically placed into the bone and creates a stable connection point capable of bearing elastic forces. It exerts its effect through mechanical fixation and load-bearing.

Relation to Medical Device Definition: IMF screws are temporary, passive implantable medical devices that operate on the principle of mechanical fixation.

➤ Tent Screw Effect Type:

The screw is mechanically fixed into the bone and creates structural support over the graft area. It exerts its effect through mechanical support and area protection.

Relation to Medical Device Definition: Tent screws are temporary, implantable passive medical devices that work on the principle of mechanical support.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

➤ Bone Fixation Screw Effect Type:

The screw mechanically fixes the graft material to the bone surface. It exerts its effect through mechanical fixation and stabilization.

Relation to Medical Device Definition: These products are passive medical devices that have no pharmacological effect and operate on the principle of mechanical fixation.

➤ Membrane Fixation Screw (Bone Tack) Effect Type:

The membrane is mechanically attached to the bone surface using bone tacks, and the regeneration area is protected. It exerts its effect through mechanical fixation and area protection.

Relation to Medical Device Definition: Bone tack products are temporary implantable passive medical devices that operate on the principle of mechanical action.

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Thanks to the smooth design of the membrane pins, it is prevented from rotating the membrane under-neath while the pin is being placed. It is possible to fix the membrane with a single blow without moving it. In addition, the additional thickness in the body part allows the short pin to come back out of its place.

Mini screws are produced with SLA (sandblasted, large grit, acid-etched) surface for faster osteointegration and stronger anchoring force.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

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Our product is a Mini Screw device, designed for long-term use and used as a surgical invasive device. For this reason, the product class has been determined as IIb due to its compliance with the article specified in Annex VIII Rule 8 of the Classification Rules of the 2017/745 / EU Medical Device Regulation and detailed below.

**Medical Device Regulation 2017/745/EU Annex IX
Full Quality Assurance System
Class IIb
2017/745/EU Annex VIII, Rule 8**

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➤ GBR Bone Fixation Screw

The bone fixation screw is intended to be used during guided bone regeneration (GBR) and similar reconstructive procedures in oral and maxillofacial surgery to **mechanically fix bone grafts or bone blocks** to the recipient bone site. The device is designed to provide stability of the graft material and to prevent unwanted movement during the healing process.

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➤ Intermaxillary Fixation Screw (IMF Screw)

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- Ridge augmentation
- Extraction sites (implant preparation/placement)
- Guided bone regeneration (GBR) procedures
- Sinus lift

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Users of these products;

- Product application includes dental implants, dental and skeletal jaw deformities, and trauma procedures. Such applications should be performed by dentists who have received the necessary training in this field.
- Must have experience in screw placement.
- Must have high attention skills during the procedure.
- Have knowledge about the use of the products.

The information provided in this manual is not sufficient for the application of screw systems and related components.

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Users of these products;

- Product application includes suitable bone substrate preparation (GBR), dental and skeletal jaw deformities and trauma procedures for dental implant application. Such applications should be carried out by dentists who have received the necessary training in this field.
- Must have experience in screw placement.
- Attention should be high during the procedure.
- Should have knowledge about the use of products.
- The information given in this user manual is not sufficient to apply screw systems and related components.

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Thanks to the smooth design of the membrane pins, it is prevented from rotating the membrane underneath while the pin is being placed. It is possible to fix the membrane with a single blow without moving it. In addition, the additional thickness in the body part allows the short pin to come back out of its place. Mini screws are produced with SLA (sandblasted, large grit, acid-etched) surface for faster osteointegration and stronger anchoring force.

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

NO	PRODUCT CODE	PRODUCT NAME	INTENDED USE
1	HSS	Screwdriver Body	A hand tool to which screw holder tips are attached for manual insertion or removal of the screw into the bone.
2	LHD	Driver (Long)	It is used to insert or remove the screw into the bone.
3	SHD	Driver (Short)	It is used to insert or remove the screw into the bone.
4	LCD	Contra Angle Driver (Long)	It is used to drive or remove the screw into the bone with a motorised hand tool.
5	SCD	Contra Angle Driver (Short)	It is used to drive or remove the screw into the bone with a motorised hand tool.
6	SM002	Drill \varnothing 1.1 mm	It is used in hard bones to create a pioneer bed for the mini screw to advance comfortably.
NO	PRODUCT CODE	PRODUCT NAME	INTENDED USE
1	HSS	Screwdriver Handle (All in One & Semi-Oss)	A hand tool to which screw holder tips are attached for manual insertion or removal of the screw into the bone.
2	LHD	Driver Tip (Long), Semi-Oss	It is used to insert or remove the screw into the bone.
3	SHD	Driver Tip (Short), Semi-Oss	It is used to insert or remove the screw into the bone.
4	LCD	Driver Tip (Contra-angle, Long), Semi-Oss	It is used to drive or remove the screw into the bone with a motorised hand tool.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

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
After the implementation and verification of all risk control measures by our company, our company decided whether the total residual risks of the medical device were acceptable in the light of the use of the criteria defined in the risk management plan.

However, within the scope of the data obtained as a result of the clinical evaluation, side effects and complications are defined in the instructions for use in line with equivalent companies and literatures.

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- ✓ Our product is disposable and should not be reused. Re-use of our product may cause infection of the patient.
- ✓ During the treatment, a maximum insertion torque of 3 Nm should be applied to the product as specified in this manual. Otherwise, the tip may break and implantation may not be realised. In addition, the patient may be harmed due to this situation.
- ✓ The shelf life of our product is determined as 5 years and the labels on the products we produce contain the production and expiry date information of our product. The patient may get an infection as a result of using the expired product.
- ✓ During the operation, the operation must be performed in accordance with the "Surgical Technique" section in this manual. In cases where these warnings are not followed, the screw may be implanted in such a way that it coincides with the patient's tooth root and as a result, implantation may not be realised. In addition, when this situation occurs, the patient may be infected and damaged.
- ✓ Surgical complications may occur during the operation. As a result of these complications, the patient may be harmed.
- ✓ Product complications may occur during the operation. As a result of these complications, treatment may not be realised.
- ✓ Side effects may occur due to the product. The patient may be harmed as a result of side effects.

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In the risk analysis study, since there is no measure that can be taken other than the usage information to be given to the user, the risk that cannot be reduced to acceptable levels;

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Damage	Hazard	Potential Cause	Reason for Non-Deductibility
Orthodontic Mini Screws			
Infection of the Patient	Biological Areas Functionality Labelling Warnings Failure to Use	Possible reuse of single-use products	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings Regarding Residual Risks" as "Our product is for single use and should not be reused. Re-use of our product may cause the patient to become infected."
Implantation does not take place Damage to the Patient	Failure to Use Operational Instructions	Tip breakage due to excessive lateral pressure	The identified risk could not be reduced to the minimum level. Since there is no precaution to be taken, it will be considered as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings on Residual Risks" under the heading "During treatment, a maximum insertion torque of 3 Nm should be applied to the product as specified in this user manual. Otherwise, the tip may break and implantation may not occur. In addition, the patient may be harmed due to the occurrence of this situation."

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/--

Infection of the Patient	Labelling Failure to Use Product Realisation Process	Use of expired products	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described under the heading "Warnings Regarding Residual Risks" of the User Manual prepared by us as "The shelf life of our product is determined as 5 years and the labels on the products we produce contain the production and expiry date information of our product. The patient may become infected as a result of using the expired product."
Implantation Failure	Failure to Use	Implanting the screw into the root of the patient's tooth	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described under the heading "Warnings Regarding Residual Risks" in the User Manual prepared by us under the heading "During the operation, the operation must be performed in accordance with the "Surgical Technique" section in this user manual. In cases where these warnings are not followed, the screw may be implanted in such a way that it coincides with the patient's tooth root and as a result, implantation may not occur. In addition, when this situation occurs, the patient may be infected and damaged."

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP
SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

Damage to the Patient	Failure to Use	Infection of the patient due to the implantation of the screw into the root of the patient's tooth	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described under the heading "Warnings Regarding Residual Risks" in the User Manual prepared by us under the heading "During the operation, the operation must be performed in accordance with the "Surgical Technique" section in this user manual. In cases where these warnings are not followed, the screw may be implanted in such a way that it coincides with the patient's tooth root and as a result, implantation may not occur. In addition, when this situation occurs, the patient may be infected and damaged." is notified to the users with the sentence.
Damage to the Patient	Operational Instructions Labelling	The occurrence of surgical complications during the operation	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings about Residual Risks" as "Surgical complications may occur during the operation. As a result of the occurrence of these complications, the patient may be harmed."
Failure to Treat	Operational Instructions Labelling	Realisation of product complications during the operation	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

			as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings on Residual Risks" as "Product complications may occur during the operation. As a result of the occurrence of these complications, treatment may not be realised."
Damage to the Patient	Warnings	Occurrence of Side Effects	The identified risk could not be reduced to the minimum level. Since there is no precaution to be taken, it will be considered as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings Regarding Residual Risks" as "Side effects may occur due to the product. The patient may be harmed as a result of the occurrence of side effects."

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The risks were evaluated as residual risks by us and a benefit/risk analysis comparison was made for residual risks as for all risks. With the comparison made, it has been concluded that our Product is usable despite the residual risks

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- ✓ The product is manufactured from pure titanium (Grade 4) raw material in accordance with EN 1642, EN 1642, TS EN ISO 19023, ISO 5832-2 and ASTM F67-13 Standard.
- ✓ Since our products are sterilised with gamma in accordance with the parameters specified in our user manuals by the users before use, the products are packed untouched in Class 100,000 Clean Room in order to ensure microbiological safety until they are sterilised.
- ✓ The gamma sterilisation method declared by us has been validated in accordance with EN 11137- 1/2 Standard.
- ✓ All special processes from raw materials to the final product have been validated by us.
- ✓ The biocompatibility of our finished products has been carried out in impartial laboratories in accordance with EN 10993-1 Standard.
- ✓ Tests for the mechanical strength of our products have been carried out in impartial laboratories in accordance with TS EN ISO 19023 Standard.

b. For Mini Screw System For Other (GBR) Purpose;

- ✓ The product is manufactured from pure titanium (Grade 4) raw material in accordance with EN 1642, EN 1642, TS EN ISO 19023, ISO 5832-2 and ASTM F67-13 Standard.
- ✓ Our products should be sterilised with autoclave before use by the users in accordance with the parameters specified in our user manuals.
- ✓ The autoclave sterilisation method declared by us has been validated in accordance with EN 17665- 1.
- ✓ All special processes from raw materials to the final product have been validated by us.
- ✓ The biocompatibility of our finished products has been carried out in impartial laboratories in accordance with EN 10993-1 Standard.
- ✓ Tests for the mechanical strength of our products have been carried out in impartial laboratories in accordance with TS EN ISO 19023 Standard.

c. For Intermxillary Fixation (IMF) Screws;

- ✓ All IMF screws are self-drilling. The risk of damaging dental roots is thus reduced as there is no need for pre-drilling.
- ✓ The product is manufactured from pure titanium (Grade 4) raw material in accordance with EN 1642, EN 1642, TS EN ISO 19023, ISO 5832-2 and ASTM F67-13 Standard.
- ✓ Our products should be sterilised with autoclave before use by the users in accordance with the parameters specified in our user manuals.
- ✓ The autoclave sterilisation method declared by us has been validated in accordance with EN 17665- 1.
- ✓ All special processes from raw materials to the final product have been validated by us.
- ✓ The biocompatibility of our finished products has been carried out in impartial laboratories in accordance with EN 10993-1 Standard.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

- ✓ Tests for the mechanical strength of our products have been carried out in impartial laboratories in accordance with TS EN ISO 19023 Standard.

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For the Dental “Mini Screw” products, equivalent devices have been identified as mini implant systems that have been widely used in clinical practice and whose performance has been well documented in the literature. The comparison demonstrated that there are no significant differences between our device and equivalent products in terms of technical, biological, and clinical aspects.

From a technical perspective, both devices are manufactured from titanium alloy with similar production methods, mechanical strength properties, and the ability to maintain long-term stability under functional loading. The screw design, implant geometry, and surface characteristics are designed at an equivalent level to support reliable stability.

From a biological perspective, both our device and equivalent mini implant systems are composed of inert and biocompatible materials suitable for long-term contact with bone and soft tissue. Neither device contains human, animal, or blood-derived components, and neither exhibits release or absorption characteristics. Therefore, there are no differences in terms of biological profile.

From a clinical perspective, mini implant systems are considered a “Well-Established Technology (WET)”. Long-term clinical evidence in the literature demonstrates their safety and effectiveness in applications such as prosthetic stabilization, orthodontic anchorage, and specific prosthetic rehabilitations. Available studies confirm that our device achieves comparable success rates, complication profiles, and patient satisfaction to equivalent products. Furthermore, post-market surveillance (PMS) data also support this equivalence.

In conclusion, our device is considered equivalent to other mini implant systems in terms of technical, biological, and clinical parameters, and is clinically acceptable regarding safety and performance.

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Our product is considered "Well-Established Technology (WET)" because it has been widely used in clinical applications for many years and its safety and performance have been extensively documented. Due to the long history of use and extensive available clinical data, the device is supported by literature, post-market surveillance (PMS) activities, and equivalence to recognized implant systems such as Straumann®. Given the well-established nature of the technology, **no specific pre-CE marking clinical investigations were conducted for this device**. Instead, the clinical evidence is derived from long-term market experience, scientific literature, and equivalent device data, all of which confirm the safety, biocompatibility, and performance of the device.

As a result, no clinical investigations prior to CE marking were required, and the available clinical evidence from WET classification, PMS data, and published literature provides sufficient demonstration of the device's clinical safety and performance.

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Clinical evidence regarding the safety and performance of the device has been obtained not only from the scientific literature but also from Post-Market Surveillance (PMS) activities, post-market clinical follow-up (PMCF) reports, user feedback, and the long-term clinical experience with similar devices. These sources demonstrate that the device has been used safely and effectively since its placement on the market. In addition, the extensive long-term clinical success of the equivalent device, Straumann® dental mini screws, as documented in the literature, provides supplementary evidence supporting the safety and performance profile of the device.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

In conclusion, the clinical data obtained from other sources, including PMS and PMCF data as well as equivalent device literature, further confirm the safety and clinical performance of the device.

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Dental **Mini Screw (Mini Vida)** systems are small-diameter implant solutions that have been widely used for many years, particularly in clinical applications such as prosthetic stabilization and orthodontic anchorage. The clinical evaluation of our product has been based on available literature, comparisons with equivalent mini implant systems, and post-market surveillance (PMS) activities. These data demonstrate that the device performs safely and effectively within its intended use and provides the expected clinical benefits in dental practice.

From a technical and biological perspective, the device is manufactured from titanium alloy, offering an inert and biocompatible structure. Its small diameter and specific design allow for minimally invasive procedures and facilitate rapid stabilization in clinical use. Reported clinical findings show that the device delivers reliable performance under functional loading, contributes to patient satisfaction, and maintains complication rates within an acceptable range.

Post-market surveillance activities have confirmed the favorable safety profile of the device, with adverse events reported at a low frequency. Supported by literature and equivalent device data, these results indicate that the device can be considered a **Well-Established Technology (WET)**.

In conclusion, the clinical performance and safety evaluation of **Mini Vida products** confirms compliance with the General Safety and Performance Requirements (GSPR) defined under the MDR. The device has been shown to ensure functional improvement, enhance patient quality of life, and provide long-term reliability, while residual risks remain within an acceptable level.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

5.5. Devam Eden Veya Planlanan Piyasaya Arz Sonrası Klinik Takip / Ongoing or Planned Post-Marketing Clinical Follow-up

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Although our products have been on the market for many years, PMCF activities are being carried out to ensure that the safety and performance of the device are continuously verified throughout its expected lifetime. Within this scope, the PMCF Plan aims to confirm the clinical benefits of the device, identify new or unforeseen side effects, monitor known side effects and contraindications, assess new risks based on real clinical evidence, and ensure the ongoing balance of benefits and risks.

The planned activities include systematic literature reviews, analysis of user feedback, evaluation of complaints and recall data, as well as continuous monitoring of clinical experience gathered from the field. PMCF activities also aim to confirm the device's conformity with its intended purpose and to identify systematic misuse or off-label use.

Based on these activities, the PMCF Report has presented up-to-date findings on the safety and performance of the device. The results demonstrated that the device achieved the expected clinical success in use, that adverse events occurred at a low rate, and that the risks remained within an acceptable level.

In conclusion, the ongoing and planned PMCF activities aim to continuously verify the safety and performance of the device, identify new risks, and demonstrate full compliance with the General Safety and Performance Requirements (GSPR) as defined under the MDR.

6. OLASI TANI VEYA TEDAVİ ALTERNATİFLERİ / POSSIBLE DIAGNOSIS OR TREATMENT ALTERNATIVES

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There are several alternative treatment methods and device options to mini screw systems. In the treatment of tooth loss, **conventional larger-diameter dental implants** are widely used in patients with sufficient bone volume. In addition, **removable full or partial dentures** are a classical option, although they are often associated with stability issues and reduced patient comfort. In orthodontic applications, alternatives to mini screws include **extraoral anchorage devices**, removable appliances, or more traditional orthodontic treatment methods. However, these approaches usually require more invasive procedures, longer treatment times, or additional appliances that may negatively impact patient compliance.

In terms of materials and surface coatings, alternatives are also available. While mini screws are most commonly manufactured from **titanium or titanium alloys**, some manufacturers prefer alternative materials such as **titanium-zirconium alloys**. Furthermore, different surface modification techniques may be applied, including **sandblasting, acid etching, plasma-sprayed coatings, or hydroxyapatite coatings**. These variations are intended to enhance biocompatibility, improve mechanical fixation with the bone, or shorten the healing period.

In conclusion, although alternative treatment methods and material/coating differences exist, mini screws are considered a reliable option due to their **minimally invasive application, ability to achieve rapid stabilization, and high clinical success rates**. They are especially preferred in cases requiring orthodontic anchorage and prosthetic stabilization.

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No specific training is required to use these products , but users of these products;

- **Product application includes suitable bone substrate preparation (GBR), dental and skeletal jaw deformities and trauma procedures for dental implant application. Such applications should be carried out by dentists who have received the necessary training in this field.**
- **Must have experience in screw placement.**
- **Attention should be high during the procedure.**
- **Should have knowledge about the use of products.**

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

➤ The information given in this user manual is not sufficient to apply screw systems and related components.

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Bkz. MDR.TD02.HLS.01 Harmonize Standart Listesi

For Orthodontic Mini Screw products manufactured and sold by us, standards, regulations, common specifications, guidelines and guidance documents accepted in our country or worldwide are taken into consideration and applied. The updates of documents such as standards, regulations, guidance documents applied for our products and system are constantly monitored and when any update occurs in the documents, this situation is examined by us and the requirements of the current documents are fulfilled.

Look: MDR.TD02.HLS.01 Harmonised Standard List

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A summary of the safety and clinical performance of the device for patients is given below.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ/SSCP SUMMARY OF SAFETY AND CLINICAL PERFORMANCE/SSCP

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Membrane Fixation Screws have a flexible driver with a flexible tip that carries them to the operation area and is used to fix them to the bone. On a hard surface, this driver is brought on the pin and pressed at an angle of 90 degrees and the pin is removed from the pin slot with a "click" sound. The screw (pin) at

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

the end of the driver is brought to a certain point of the membrane to be fixed and sent into the bone by hand force or by hitting the lower part of the driver with a hammer. After the screw is inserted into the jawbone, the driver is slightly rotated to the right and left and the screw is left on the membrane by stretching the driver tip. After the time required for GBR (3-10 months) is completed, these screws (pins) are usually removed. Since they are made of biocompatible material, some physicians do not remove these pins at all and leave them in place. There is no medical disadvantage in leaving the pins in the mouth.

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Bone Fixation Mini Screws are used to fix the block bone graft to the donor site by drilling beforehand. Both the block bone graft and the bone in the donor site are drilled simultaneously with the drill suitable for the screw diameter. The bone fixation screw is inserted into the prepared hole by turning clockwise with the help of a driver (screwdriver) and the block bone graft is stabilized in the desired area. Depending on the size of the bone graft, more than one screw can be used to prevent the graft from moving. There are driver tips that can be used manually and with contra-angle handpiece and are the same as those used in tent screws. If the contra-angle handpiece compatible tip is preferred as the driver during placement, the screws are sent clockwise (forward) into the bone using a speed of 30 rpm and a torque of 10 Ncm.

Intermaxillary Fixation (IMF) Screws. All IMF screws are self-drilling. The risk of damaging dental roots is thus reduced as there is no need for pre-drilling. The ideal position is above the teeth roots. Using the screwdriver, the screw is directly inserted into the bone; pre-drilling is not necessary. The screw is advanced until the concave shoulder of the screw head is 1 mm from the mucosal surface. The screw may be readjusted for proper alignment of the wire holes contained within the screw head. As an alternative to the above procedure, the mini screw may be also inserted by using the contra-angle handpiece. In case of insertion by using a contra-angle handpiece, a speed of 30 rpm and a torque of 30 Ncm must be set.

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➤ **GBR Bone Fixation Screw**

The bone fixation screw is intended to be used during guided bone regeneration (GBR) and similar reconstructive procedures in oral and maxillofacial surgery to **mechanically fix bone grafts or bone blocks** to the recipient bone site.

The device is designed to provide stability of the graft material and to prevent unwanted movement during the healing process.

➤ **GBR Tenting Screw**

The tenting screw is intended for use in guided bone regeneration (GBR) procedures to prevent collapse of the soft tissue over the bone graft and to **maintain the space required for bone regeneration**. The device is placed into the bone surface to support the membrane or soft tissue and to create an adequate space for new bone formation.

➤ **Membrane Fixation Screw (Bone Tack)**

The membrane fixation screw (bone tack) is intended to be used to **fix and stabilize barrier membranes** applied during guided bone regeneration (GBR) and guided tissue regeneration (GTR) procedures onto the recipient bone surface.

The device prevents undesired movement of the membrane at the surgical site and enhances the predictability of the regenerative process.

➤ **Orthodontic Mini Screw**

The orthodontic mini screw is intended to provide **temporary skeletal anchorage** during orthodontic treatments in order to facilitate controlled tooth movement.

The device is inserted into the alveolar bone to support the application of orthodontic forces and is intended to be removed after completion of the treatment.

➤ **Intermaxillary Fixation Screw (IMF Screw)**

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

The intermaxillary fixation screw (IMF screw) is intended to be used for the **temporary fixation of the upper and lower jaws** during the treatment of jaw fractures, orthognathic surgery, and post-traumatic conditions.

The device allows intermaxillary fixation through the use of elastics or wire systems.

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Block Bone Graft Fixation Screw

- Stabilization and fixation of block bone grafts in oral and maxillofacial surgery.
- Used for bone grafts in the oral cavity, providing stability for guided bone regeneration (GBR).
- Indicated in cases requiring ridge augmentation, alveolar reconstruction, and graft fixation.
- Supports intermaxillary fixation when used adjunctively in certain surgical procedures.

Tent Screw (Tent Peg Screw)

- The Tenting Screw is used to stabilize, fixate, and/or support bone grafts, bone filling materials and /or barrier membranes used for regeneration of bone in the oral cavity

Membrane Fixation Screw / Pin

- Periodontal / intrabony defects
- Ridge augmentation
- Extraction sites (implant preparation/placement)
- Guided bone regeneration (GBR) procedures
- Sinus lift

Orthodontic Mini Screw (Temporary Anchorage Device – TAD)

- Provides temporary anchorage in orthodontic treatments.
- Indicated for:
 - Treatment of Class I, Class II, and Class III malocclusion.
 - Anchor recovery and anchor reinforcement.
 - Asymmetric case management.
 - Steepening of upper and lower molars.
 - Correction of overbite (molars, premolars, incisors).
 - Management of deep bite and open bite conditions.
- Reduces reliance on patient compliance for orthodontic mechanics.

Intermaxillary Fixation Screw (IMF Screw)

- Indicated for intermaxillary fixation in oral and maxillofacial surgery.
- Used in fracture management of the mandible and maxilla.
- Provides temporary stabilization of occlusion during healing.
Can also assist in bone graft stabilization in combined surgical cases.

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Users of these products;

- Product application includes dental implants, dental and skeletal jaw deformities, and trauma procedures. Such applications should be performed by dentists who have received the necessary training in this field.
 - Must have experience in screw placement.
 - Must have high attention skills during the procedure.
 - Have knowledge about the use of the products.
- The information provided in this manual is not sufficient for the application of screw systems and related components.

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- Patients with no or inadequate bone support at the insertion site.
- Foreign material sensitisation. Where material sensitivity is suspected, appropriate testing should be performed prior to implantation to determine if sensitivity is present.
- Conditions that tend to impair the patient's ability or willingness to cooperate and follow instructions during the recovery period.
- Degenerative diseases whose progression may adversely affect implant placement.
- Insufficient amount of bone for secure implant anchorage.
- Blood flow restrictions, radiation therapy and/or previous infections that may delay healing and increase the likelihood of infection and/or rejection of the implant.
- Inadequate coverage with healthy tissue.
- Procedures where there is a non-sterile environment, i.e. open cavities such as sinuses.

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Thanks to the thin diameters of 1 mm and 1.2 mm of bone fixation screws, breakage or cracking of the block graft is prevented. However, low volume formation in the donor area positively affects the success of the GBR (guided bone regeneration) method, resulting in greater volume of bone formation.

The most important advantage of tenting screws over conventional screws is that, thanks to their large head, they act as a tent under the membrane and prevent the membrane from collapsing into the graft particles. In addition, possible secondary infections are prevented thanks to the shiny surface on the body part.

Thanks to the smooth design of the membrane pins, it is prevented from rotating the membrane under-neath while the pin is being placed. It is possible to fix the membrane with a single blow without moving it. In addition, the additional thickness in the body part allows the short pin to come back out of its place. Mini screws are produced with SLA (sandblasted, large grit, acid-etched) surface for faster osteointegration and stronger anchoring force.

The head of the intermaxillary fixation mini screws is produced in a large and convex shape to prevent soft tissue trauma such as lip-cheek during the temporary period in the mouth.

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Previous variants of the device and are not available

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All our Mini Screw models can be used with the following hand tools.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

NO	PRODUCT CODE	PRODUCT NAME	INTENDED USE
1	HSS	Screwdriver Body	A hand tool to which screw holder tips are attached for manual insertion or removal of the screw into the bone.
2	LHD	Driver (Long)	It is used to insert or remove the screw into the bone.
3	SHD	Driver (Short)	It is used to insert or remove the screw into the bone.
4	LCD	Contra Angle Driver (Long)	It is used to drive or remove the screw into the bone with a motorised hand tool.
5	SCD	Contra Angle Driver (Short)	It is used to drive or remove the screw into the bone with a motorised hand tool.
6	SM002	Drill \varnothing 1.1 mm	It is used in hard bones to create a pioneer bed for the mini screw to advance comfortably.
NO	PRODUCT CODE	PRODUCT NAME	INTENDED USE
1	HSS	Screwdriver Handle (All in One & SemiOss)	A hand tool to which screw holder tips are attached for manual insertion or removal of the screw into the bone.
2	LHD	Driver Tip (Long), Semi-Oss	It is used to insert or remove the screw into the bone.
3	SHD	Driver Tip (Short), Semi-Oss	It is used to insert or remove the screw into the bone.
4	LCD	Driver Tip (Contra-angle, Long), Semi-Oss	It is used to drive or remove the screw into the bone with a motorised hand tool.



GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

5	SCD	Driver Tip (Contra-angle, Short), Semi-Oss	It is used to drive or remove the screw into the bone with a motorised hand tool.
6	SM002	Drill Ø1.1 mm (Semi-Oss)	It is used in hard bones to create a pioneer bed for the mini screw to advance comfortably.

12.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Do-
cument No

RP.18.03

QMS Template Release
Date


16.08.2021

QMS Template
Revision No

00

QMS Template Revision
Date

--/--/---

13.

14.

NO	PRODUCT CODE	PRODUCT NAME	INTENDED USE
1	HSS	Screwdriver Body	A hand tool to which screw holder tips are attached for manual insertion or removal of the screw into the bone.
2	SMHD	GBR Kit Driver (Long)	It is used to insert or remove the screw into the bone.
3	SMT009	GBR Kit Driver (Short)	It is used to insert or remove the screw into the bone.
4	SMT010	GBR Kit Contra Angle Driver	It is used to drive or remove the screw into the bone with a motorised hand tool.
5	SMT006	Pin Holder (Straight)	Specially produced for manual insertion or removal of the screw into the bone, it acts as a holder on the end.
6	SMT007	Pin Holder (Angled)	Specially produced for manual insertion or removal of the screw into the bone, it acts as a holder on the end.
7	SMT011	GBR Kit Drill	It is used in hard bones to create a pioneer bed for the mini screw to advance comfortably.

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
GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP

SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/--

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

NO	PRODUCT CODE	PRODUCT NAME	INTENDED USE
1	HSS	Screwdriver Handle (All in One & Semioss)	A hand tool to which screw holder tips are attached for manual insertion or removal of the screw into the bone.
2	SMHD	Driver Tip (Long) All in One	It is used to insert or remove the screw into the bone.
3	SMT009	Driver Tip (Short) All in One	It is used to insert or remove the screw into the bone.
4	SMT010	Driver Tip (Contra-angle), All in One	It is used to drive or remove the screw into the bone with a motorised hand tool.
5	SMT006	Pin Holder (Straight) (All in One & MPK)	Specially produced for manual insertion or removal of the screw into the bone, it acts as a holder on the end.
6	SMT007	Pin Holder (Angled) (All in One & MPK)	Specially produced for manual insertion or removal of the screw into the bone, it acts as a holder on the end.
7	SMT011	Drill Ø1 mm (All in One)	It is used in hard bones to create a pioneer bed for the mini screw to advance comfortably.

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Residual / Residual Risk: Risks that cannot be reduced to the level of "acceptable risk" despite all measures taken against the risks that occur or may occur related to the products are defined as residual risks. Below are warnings about the residual risks of our product.

- ✓ Our product is disposable and should not be reused. Re-use of our product may cause infection of the patient.
- ✓ During the treatment, a maximum insertion torque of 3 Nm should be applied to the product as specified in this manual. Otherwise, the tip may break and implantation may not be realised. In addition, the patient may be harmed due to this situation.
- ✓ The shelf life of our product is determined as 5 years and the labels on the products we produce contain the production and expiry date information of our product. The patient may get an infection as a result of using the expired product.
- ✓ During the operation, the operation must be performed in accordance with the "Surgical Technique" section in this manual. In cases where these warnings are not followed, the screw may be implanted in such a way that it coincides with the patient's tooth root and as a result, implantation may not be realised. In addition, when this situation occurs, the patient may be infected and damaged.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/--

- ✓ Surgical complications may occur during the operation. As a result of these complications, the patient may be harmed.
- ✓ Product complications may occur during the operation. As a result of these complications, treatment may not be realised.
- ✓ Side effects may occur due to the product. The patient may be harmed as a result of side effects.

Proimtech Sağlık Ürünleri A.Ş., as the manufacturer and seller of these products, is not liable for any direct or indirect damages resulting from misuse or improper use or application, in particular as a result of improperly directed and unsuitable operations and service. If unauthorized persons undertake repair, replacement, or use of Bilimplant products in conjunction with products of other manufacturers, the manufacturer cannot be held responsible for the consequences; furthermore, warranty rights are reserved.

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Damage	Hazard	Potential Cause	Reason for Non-Deductibility
Orthodontic Mini Screws			
Infection of the Patient	Biological Areas Functionality Labelling Warnings Failure to Use	Possible reuse of single-use products	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings Regarding Residual Risks" as "Our product is for single use and should not be reused. Reuse of our product may cause the patient to become infected."

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP
SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

<p>Implantation does not take place Damage to the Patient</p>	<p>Failure to Use Operational Instructions</p>	<p>Tip breakage due to excessive lateral pressure</p>	<p>The identified risk could not be reduced to the minimum level. Since there is no precaution to be taken, it will be considered as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings on Residual Risks" under the heading "During treatment, a maximum insertion torque of 3 Nm should be applied to the product as specified in this user manual. Otherwise, the tip may break and implantation may not occur. In addition, the patient may be harmed due to the occurrence of this situation."</p>
<p>Infection of the Patient</p>	<p>Labelling Failure to Use Product Realisation Process</p>	<p>Use of expired products</p>	<p>The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described under the heading "Warnings Regarding Residual Risks" of the User Manual prepared by us as "The shelf life of our product is determined as 5 years and the labels on the products we produce contain the production and expiry date information of our product. The patient may become infected as a result of using the expired product."</p>

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/--

Implantation Failure	Failure to Use	Implanting the screw into the root of the patient's tooth	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described under the heading "Warnings Regarding Residual Risks" in the User Manual prepared by us under the heading "During the operation, the operation must be performed in accordance with the "Surgical Technique" section in this user manual. In cases where these warnings are not followed, the screw may be implanted in such a way that it coincides with the patient's tooth root and as a result, implantation may not occur. In addition, when this situation occurs, the patient may be infected and damaged."
Damage to the Patient	Failure to Use	Infection of the patient due to the implantation of the screw into the root of the patient's tooth	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described under the heading "Warnings Regarding Residual Risks" in the User Manual prepared by us under the heading "During the operation, the operation must be performed in accordance with the "Surgical Technique" section in this user manual. In cases where these warnings are not followed, the screw may be implanted in such a way that it coincides with the patient's tooth root and as a result, implantation may not occur."

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

			<p>cur. In addition, when this situation occurs, the patient may be infected and damaged." is notified to the users with the sentence.</p>
Damage to the Patient	Operational Instructions Labelling	The occurrence of surgical complications during the operation	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings about Residual Risks" as "Surgical complications may occur during the operation. As a result of the occurrence of these complications, the patient may be harmed."
Failure to Treat	Operational Instructions Labelling	Realisation of product complications during the operation	The identified risk could not be reduced to the minimum level. Since there are no measures to be taken, it will be considered as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings on Residual Risks" as "Product complications may occur during the operation. As a result of the occurrence of these complications, treatment may not be realised."
Damage to the Patient	Warnings	Occurrence of Side Effects	The identified risk could not be reduced to the minimum level. Since there is no precaution to be taken, it will be considered

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

as a residual risk by us. This residual risk is described in the User Manual prepared by us under the heading "Warnings Regarding Residual Risks" as "Side effects may occur due to the product. The patient may be harmed as a result of the occurrence of side effects."

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- ✓ All IMF screws are self-drilling. The risk of damaging dental roots is thus reduced as there is no need for pre-drilling.
- ✓ The product is manufactured from pure titanium (Grade 4) raw material in accordance with EN 1642, EN 1642, TS EN ISO 19023, ISO 5832-2 and ASTM F67-13 Standard.
- ✓ Our products should be sterilised with autoclave before use by the users in accordance with the parameters specified in our user manuals.
- ✓ The autoclave sterilisation method declared by us has been validated in accordance with EN 17665- 1.
- ✓ All special processes from raw materials to the final product have been validated by us.
- ✓ The biocompatibility of our finished products has been carried out in impartial laboratories in accordance with EN 10993-1 Standard.
- ✓ Tests for the mechanical strength of our products have been carried out in impartial laboratories in accordance with TS EN ISO 19023 Standard.

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From a clinical perspective, mini implant systems are considered a “**Well-Established Technology (WET)**”. Long-term clinical evidence in the literature demonstrates their safety and effectiveness in applications such as prosthetic stabilization, orthodontic anchorage, and specific prosthetic rehabilitations. Available studies confirm that our device achieves comparable success rates, complication profiles, and patient satisfaction to equivalent products. Furthermore, post-market surveillance (PMS) data also support this equivalence.

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Our product is considered “**Well-Established Technology (WET)**” because it has been widely used in clinical applications for many years and its safety and performance have been extensively documented. Due to the long history of use and extensive available clinical data, the device is supported by literature, post-market surveillance (PMS) activities, and equivalence to recognized implant systems such as Straumann®.

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

Given the well-established nature of the technology, **no specific pre-CE marking clinical investigations were conducted for this device.** Instead, the clinical evidence is derived from long-term market experience, scientific literature, and equivalent device data, all of which confirm the safety, biocompatibility, and performance of the device.

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Clinical evidence regarding the safety and performance of the device has been obtained not only from the scientific literature but also from Post-Market Surveillance (PMS) activities, post-market clinical follow-up (PMCF) reports, user feedback, and the long-term clinical experience with similar devices. These sources demonstrate that the device has been used safely and effectively since its placement on the market.

In addition, the extensive long-term clinical success of the equivalent device, Straumann® dental mini screws, as documented in the literature, provides supplementary evidence supporting the safety and performance profile of the device.

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The planned activities include systematic literature reviews, analysis of user feedback, evaluation of complaints and recall data, as well as continuous monitoring of clinical experience gathered from the field. PMCF activities also aim to confirm the device's conformity with its intended purpose and to identify systematic misuse or off-label use. Based on these activities, the PMCF Report has presented up-to-date findings on the safety and performance of the device. The results demonstrated that the device achieved the expected clinical success in use, that adverse events occurred at a low rate, and that the risks remained within an acceptable level.

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- **Product application includes suitable bone substrate preparation (GBR), dental and skeletal jaw deformities and trauma procedures for dental implant application. Such applications should be carried out by dentists who have received the necessary training in this field.**
- **Must have experience in screw placement.**
- **Attention should be high during the procedure.**
- **Should have knowledge about the use of products.**
- **The information given in this user manual is not sufficient to apply screw systems and related components.**

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02	05.01.2026	Document has been revised in light of detailed technical documentation review conducted by the notified body.	<input type="checkbox"/> Evet /Yes Geçerli Kılınan Dil /Validated Language: <input type="checkbox"/> Hayır / No
----	------------	---	---

Hazırlayan Yapan/Onay
Prepared by /Approval

Kontrol Eden/Onay
Controlled by /Approval

GÜVENLİK VE KLİNİK PERFORMANS ÖZETİ RAPORU-SSCP SAFETY AND CLINICAL PERFORMANCE SUMMARY REPORT-SSCP

QMS Template Document No	QMS Template Release Date	QMS Template Revision No	QMS Template Revision Date
RP.18.03	16.08.2021	00	--/--/---

GÜVENLİLİK VE KLİNİK PERFORMANS ÖZETİ/SSCP SUMMARY OF SAFETY AND CLINICAL PERFORMANCE/SSCP

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This summary of safety and clinical performance (SSCP) is intended to provide public access to an updated summary of key aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more comprehensive summary of the safety and clinical performance of the device for healthcare professionals is included in the first section of this document.

The SSCP is not intended to provide general advice on the treatment of a medical condition. If you have any questions about your medical condition or the use of the device in your condition, please contact your healthcare professional. This SSCP is not intended to replace an implant card or instructions for use to provide information on the safe use of the device.